

## DATA QUALITY SHEET



### **Anti-Taq monoclonal Antibody**

FOR RESEARCH USE ONLY

Cat. GC-029-0100, GC-029-0250, GC-029-0500, GC-029-1000

<b>DESCRIPTION</b>	<p>The Taq Start Monoclonal Antibodies were derived from a hybridoma (fusion of mouse myeloma cell and the cells after mouse immunization with <i>Taq</i> DNA Polymerase). Hot Start <i>Taq</i> Monoclonal Antibody is mouse IgG 2b isotype.</p> <p>Hot Start <i>Taq</i> Monoclonal Antibody inhibits polymerase activity before the onset of thermal cycling, preventing nonspecific amplification and primer-dimer formation. When the temperature is raised, the antibody is quickly inactivated and PCR proceeds. Hot Start <i>Taq</i> Monoclonal Antibody is effective with a variety of commercially available <i>Taq</i> DNA polymerases (native or recombinant). The use of Hot Start <i>Taq</i> Monoclonal Antibody significantly improves the specificity of PCR amplification what is especially important for PCR-based diagnostics.</p>
<b>APPLICATION</b>	<ul style="list-style-type: none"><li>- Hot Start PCR</li><li>- PCR-based diagnostics</li></ul>
<b>CONCENTRATION</b>	4 mg/ml
<b>UNIT DEFINITION</b>	One unit is defined as the amount of Hot Start Taq Monoclonal Antibody required to block 50% of 1 µg Taq polymerase at 37°C.
<b>U/mg RATIO</b>	2500 units of specific activity are equal to 1 mg of antibodies
<b>STORAGE BUFFER</b>	10 mM Tris-HCl (pH 7.0 at 22°C), 50 mM KCl, 0.1 mM EDTA, 50% glycerol
<b>STORAGE TEMP.</b>	Store Hot Start Taq Monoclonal Antibody at -20°C in a constant temperature freezer.
<b>REACTION BUFFER</b>	The Hot Start Taq Monoclonal Antibody reaction buffer is the same buffer used for the thermostable DNA polymerase.
<b>PURITY</b>	More than 95% in SDS electrophoresis in 15% PAAG
<b>ASSOC. ACTIVITIES</b>	No conversion to the covalently closed circular DNA to the nicked or linear form was observed after incubation of 1 µg of pUC19 with antibodies in final concentration of 6 u/µl in 20 µl of reaction mixture containing 25 mM Tris-HCl (pH 7.9), 100 mM NaCl, 10 mM MgCl <sub>2</sub> after 16 hours at 37°C.
<b>PROTOCOL</b>	The amount of antibodies required for <i>Taq</i> polymerase activity inhibition depends not on the units of the enzyme, rather the amount of <i>Taq</i> polymerase as protein! The ratio units/mg of <i>Taq</i> polymerase varies strongly from preparation to preparation of the different producers. We recommend to use 2 µl of antibody per 8 µl (40 units) of BioTherm™ <i>Taq</i> DNA polymerase.
<b>SHELF LIFE</b>	2 years from date of receipt under proper storage conditions.
<b>FUNCTIONAL ANALYSIS</b>	Tested functionally in a unit activity test.